Get to know also our other products

SaguaroX

Wireless remote controlled translation stages for easy CT sample mounting and centring.



Phantoms

Phantoms for nano CT measurements (voxel size calibration), precise spatial resolution estimation (ASTM E1695-95, nano and micro CT) and other.



Accessories

CT hedgehogs, magnetic adapters, rod-holders, rigid bases, and other accessories.

Visit www.cactux.cz also for the technical parameters of our products.



About Cactux s.r.o.

Being based in Brno, the Czech Republic, CactuX is a spin-off company founded by experienced researchers from the Laboratory of X-ray micro and nano computed tomography from CEITEC Brno University of Technology.

We develop, produce, and sell innovative accessories for industrial micro and nano CTs.

Cactux s.r.o. Hlavní 104 C 664 31 Lelekovice Czech Republic



Sales and prices: sales@cactux.cz + 420 771 281 403



Essential Toolkit NANO

Complete set of all necessary tools for nano CT sample preparation

NEW



Do you want to purchase or sell this product? Contact us!



What is

Essential Toolkit NANO?

- A complete set of all necessary tools for nano CT sample preparation.
- Contains quality step rods, tailpieces, scalpel, needle, tweezers, scissors, glue, and other tools.
- Suited for working with various sample types, ranging from firm to powdery; to be used for fixating, slicing, manipulating, or holding samples, and for a flexible and strong fixation.
- Completed by a team of experienced CT-lab technicians.
- By purchasing the kit, you will acquire all the necessary tools, all at once and without the need of buying them from different sources. When buying a new nano CT, the kit also solves the issue of identifying the essential tools used for the sample preparation.
- As it is packed in a design box, it may serve as a practical and unique business gift.

All tools necessary for nano CT sample preparation in one place. Do you want to purchase or sell this product?

Get your price quote!



Scan to get your price quote.



sales@cactux.cz: + 420 771 281 403

What does the kit contain? Technical information

Standard rods	Diameter (mm)	Quantity
Cylindrical-shaped metallic holders for the most sample's types. Length = 55 mm.	1.5	1
	3.0	1
	5.0	1

Specialized rods	Tool type	Diameter (mm)	Quantity
Specially-shaped	Step-rod	3.0	1
metallic holders compatible with	Narrowed	1.0 (top)	1
Kapton tubes or allowing easy attachment of long, thin samples. Length = 55 mm.	1.5 (bottom)	'	
		2.0 (top)	,
		3.0 (bottom)	'

Tailpieces	Diameter (mm)	Quantity
	1.5	1
For holding cylindrical standard or specialized rods.	3.0	1
	5.0	1

Sample attachment	Tool type	Quantity
Tools allowing flexible or strong fixation.	Double-sided tape	1
	Superglue	1
	Dental wax	1

Sample insulation	Diameter (mm)	Quantity
Tubing tools for powdery or fragile samples preparation and their insulation from airflows in	1.0	1
measurement chamber. Length = 5 mm.	2.0	1

Other instruments	Tool type	Quantity
	Scalpel	1
	Spare blades	5
For sample's slicing and manipulation.	Tweezers	1
mamparation.	Scissors	1
	Preparation needle	1